

Measurement form F1.1

Upholstered one piece back for Classic and Tilt

Choose this kind of cushion for the Kiddo Classic or Kiddo Tilt if you are looking for a comfortable custom made back (anatomical shaped).


This one piece back is upholstered in hard-wearing fabric on the sides and a breathable fabric on the top. Other fabrics can be chosen in the order form.

The right size for the back is automatically created from the specified dimensions of the backrest height and seat width on the order form.
No need to fill in the measurements!




If you want a different kind of cushion? Then fill in one of the other measurement forms below.




Measurement form F1.2
Upholstered one piece back with thoracic side supports



Measurement form F2
Upholstered two split back with thoracic side supports



Measurement form F3
Deep contour back



Size	Description	Value (cm)
This measurement form is part of:		
H	Total height of the back	
T	Height of the tulip <i>(can be zero)</i>	
Bp	Width of the inside of the back	
Dp	Depth of the inside of the back	