


# Measurement form F10.1

## Anatomic seat for Space and Up

Choose this cushion for the Kiddo Space or Kiddo Up if you are looking for a comfortable cushion. Light anatomic formed seat cushion with light threshold (1cm). Abduction at the front and curved upholstery at the sides




This seat cushion is upholstered in black, hard-wearing fabric on the sides and a breathable fabric on the top.

**No need to fill in the measurements!** The right size for the cushion is automatically created from the specified dimensions of the seat depth and seat width on the order form.

If you want a different kind of seat cushion? Then fill in one of the other measurement forms below.




Measurement form: F10.2  
Abduction wedge / pommel



Measurement form F10.3  
Ergo-seat cushion



Measurement form: F10.4  
Saddle seat / abduction pommel



Size	Description	Value (cm)
B	Total width of the seat cushion	= Seat width (orderform)
L	Diepte van het complete zitkussen	= Seat depth + 4 cm (orderform)