

Abduction wedge for Space and Up

Choose this option for your anatomic seat or ergo seat cushion when you are looking for an abduction wedge in order to improve the position of the legs.


The abduction wedge is upholstered in breathable fabric (black).

Please fill in the table below and add the form to the order. Don't forget to also fill in the right measurement form for the requested seat cushion!



If you want a different kind of seat cushion? Then fill in one of the other measurement forms below.




Measurement form: F10.1
Anatomic seat cushion



Measurement form F10.3
Ergo-seat cushion



Measurement form: F10.4
Saddle seat / abductionommel



Size	Description	Value (cm)
B	Total width of the seat cushion	= Seat width (orderform)
L	Total length of the seat cushion	= Seat depth + 4 cm (orderform)
Bv	Width at front side of the wedge	
Ba	Width at the rear of the wedge	Normal standard determined, but if desired, fill in value.
La	Total length of the wedge	
H	Height abduction wedge above seat	
T	Length of the wedge in front of the seat cushion	