


## Upholstered two split back with thoracic side supports for Classic and Tilt

Choose this kind of back for the Kiddo Classic or Kiddo Tilt if you are looking for a comfortable custom made back with additional support for the upper body to improve the seating position by adjusting the left and right support separately.

This one piece back is upholstered in hard-wearing fabric on the sides and a breathable fabric on the top. Other fabrics can be chosen in the order form.

Please fill in the table below and add the form to the order.




If you want a different kind of cushion? Then fill in one of the other measurement forms below.



Measurement form F1.1  
Upholstered one piece back



Measurement form F1.2  
Upholstered one piece back with thoracic side supports



Measurement form F3  
Deep contour back



Size	Description	Left support	Right support
This measurement form is part of:			
H	Total height of the back side support		
B	Width of back side support		
L	Depth of back side support		
A	Difference between height back and side support (can be zero)		
Dr	Height of the side support (back) (maximum H-A)		
Db	Height of the side support (breast side) (maximum Dr-4cm)		