


Measurement form F3

Deep contour back for Classic and Tilt

Choose this kind of back for the Kiddo Classic or Kiddo Tilt if you are looking for a comfortable custom made back with additional support for the entire upper body to improve the seating position.

This one piece back is upholstered in hard-wearing fabric on the sides and a breathable fabric on the top. Other fabrics can be chosen in the order form.

Please fill in the table below and add the form to the order.




If you want a different kind of cushion? Then fill in one of the other measurement forms below.





Measurement form F1.1
Upholstered one piece back



Measurement form F1.2
Upholstered one piece back with thoracic side supports



Measurement form F2
Upholstered two split back with thoracic side supports



Size	Description	Value (cm)
This measurement form is part of:		
H	Total height of the back	
T	Height of the tulip <i>(can be zero)</i>	
Bc	Width of the inside of the back contour	
Dc	Depth of the inside of the back contour	
C	Height of the back contour	