


Measurement form F4

Standard anatomic seat for Classic and Tilt

Choose this kind of seat for the Kiddo Classic or Kiddo Tilt if you are looking for a comfortable custom made seat. Light anatomic formed seat with light threshold (1cm). Abduction at the front and curved upholstery at the sides.




This anatomic seat is upholstered in hard-wearing fabric on the sides and a breathable fabric on the top. Other fabrics can be chosen in the order form.

The right size for the seat is automatically created from the specified dimensions of the seat width and seat depth on the order form.


No need to fill in the measurements!

If you want a different kind of back? Then fill in one of the other measurement forms below.




Measurement form F5.1 / 5.2
Ergo seat or deep contour seat



Measurement form F6.1 / 6.2
Abduction wedge for Classic or Tilt



Measurement form F7
Saddle seat / abduction seat



| Size | Description | Value (cm) |
|------|--------------------------|---------------------------------|
| B | Total width of the seat | = Seat width (order form) |
| L | Total length of the seat | = Seat depth + 6cm (order form) |