


# Measurement form F5.1


## Ergo seat for Classic and Tilt

Choose this kind of seat for the Kiddo Classic or Kiddo Tilt if you are looking for a comfortable custom made seat. Prevent sliding forwards by determining the height and length of the threshold.




If you want a different kind of back? Then fill in one of the other measurement forms below.


Measurement form F4  
Anatomic seat




Measurement form F5.2  
Deep contour seat



Measurement form F6.1 / F6.2  
Abduction wedge for Classic or Tilt





Measurement form F7  
Saddle seat



This ergo seat is upholstered in hard-wearing fabric on the sides and a breathable fabric on the top. Other fabrics can be chosen in the order form.

Please fill in the table below and add the form to the order.



Size	Description	Value (cm)
This measurement form is part of:		
B	Total width of the seat	= Seat width (order form)
L	Total length of the seat	= Seat depth + 6cm (order form)
Ld	Length of the threshold	
Hd	Height of the threshold	