

Measurement form F9.2

One piece back cushion and thoracic side supports for Space and Up

Choose this kind of cushion for the Kiddo Space or Kiddo Up if you are looking for a comfortable cushion with additional support for the upper body to improve the seating position.


This One piece back cushion and thoracic side supports are upholstered in hard-wearing fabric on the sides and a breathable fabric on the top.

Please fill in the table below and add the form to the order.





If you want a different kind of cushion? Then fill in one of the other measurement forms below.

Measurement form: F9.1
Standard, one piece back cushion



Measurement form F9.3
One piece back, deep contour cushion



Size	Description	Value (cm)		
H	Total height of the back cushion	= Backrest height (order form)		
T	Height of the tulip	ZB 20 - 22	ZB 22 - 31	ZB 31 - 36
		3,5	4,5	5,5
A	Difference top of the seat and side support			
B	Width between the side supports			
L	Depth of the side support			
Dr	Height of the side support (back)			
Db	height of the side support (breast side)			