

Specifications

For the Kiddo Up the following data are applicable.

	Seat width 16 - 32 cm		Overall length from 65 cm
	Seat depth 16 - 38 cm		Overall width from 42 cm
	Seat height 34 - 50 cm		Back angle 10 degrees
	Leg length 15 - 40 cm		Camber 7 degrees
	Back height 16 - 40 cm		Max. user weight 50 kg.
	Kg. Overall weight from 9 kg included complete seat cushion and back cushion		Rear wheel size 20" 22" 24"
	Kg. Weight without rear heels from 6,4 kg included complete seat cushion and back cushion		Front wheel size 4" 5"

Standardization

Maximum user weight is 50kg. The Kiddo Up has been tested and approved by TÜV Rheinland for NEN-EN 12182 and NEN-EN 12183. Veldink4kids is affiliated with Nationaal Keurmerk Hulpmiddelen (NL)



MURRAY'S

Up
Kiddo by Veldink

The Kiddo Up is an active lightweight wheelchair with a tilt function especially for the smallest users. Extremely light in weight but retaining the familiar Veldink4kids quality. The Kiddo Up is always manufactured to the measurements of the child, custom made!

Do you want to contact us?



Block C Airton Park,
Airton Rd, Tallaght,
Dublin 24



Eire Tel: 01 855 5696
UK Tel: 028 91 877951



info@murrays.ie
sales@murraysmedical.co.uk

Do you want to follow us?







www.facebook.com/murrays.ie
www.facebook.com/murraysuk



www.twitter.com/murraysmed
www.twitter.com/murraysmeduk



www.murrays.ie
www.murraysmedical.co.uk

-  Tilt function for multiple positions
-  Lightweight wheelchair
-  For early childhood rehab
-  Perfectly adjustable to each user's therapeutic and physical demands

 **Veldink4kids**
no limits

Kiddo Up

The Kiddo Up is an active lightweight wheelchair with a tilt function especially for the smallest users. The tilt function offers the ability to alternate position. An active or passive position in one wheelchair.



An innovative wheelchair for kids who can hoop for themselves and sometimes need a tilt function!

Foldable pushbar

Fully personalize the wheelchair by the color, wheel protectors and back image designs.

One piece footplate



The Kiddo Up is always provided with a foldable one piece footplate. The one piece footplate can be provided in a heel edge version to prevent backward movement of the feet. The footplate is height and angle adjustable.

Rear wheels

Choice of 20, 22 or 24" Quick Release rear wheels with tyres of your choice. Solid or air tyres in standard or outdoor styles.

Brakes

Brakes integrated into the mudguards for easy access for the child with optional lock or attendant controlled.

Rear wheel with alu hoop




Outdoor tyre



Seats

We have a wide range of seating options to suit children's varying needs. Each type of seat or backrest cushion is individually customised and made to measure for the child with a wide range of options and additional support accessories.



Abduction wedge

Anatomic seat

Ergo seat

Pony seat

One piece back

Side support

Support shaped

Spoke protectors

Fully personalise your wheelchair with bespoke spokeguards and trays with virtually unlimited possibilities. We can even incorporate your own images or photos.



Choose your own spoke protectors! www.murrays.ie or www.murraysmedical.co.uk


Headrests

Our different types of headrests ensure optimal position of the head in combination with the right support. The headrest is adjustable in angle, height and depth. Of course we fabricate, if desired, fully customized headrests. By choosing for a stainless steel version, the strength of the position is guaranteed.



Work trays

The work trays are provided with a double bar to be mounted on the fenders. This system is easy to adjust in height and lock on the right height.



Removable work tray, depth and height adjustable. Equipped with double insertion tube for placement in the tubes on the front.

More options

There are more options possible than displayed in this flyer. Check:

www.murrays.ie
www.murraysmedical.co.uk
www.veldink4kids.com

and the comprehensive option brochures for all the possible options for the Kiddo Up.