

Accessories / Options

Legrest options
suitable for X3, X4, X5



Armrest options
suitable for X2, X3, X4, X5



Universal Headrest



All Aktiv Wheelchairs have:

- Lightweight Aluminium Frames
- Folding Frame for Easy Storage
- Crash Test Passed for Occupied Transportation
- 5 Year Frame Warranty
- Aluminium (Not Plastic) in Stress Areas
- Flip-Up Backrest, Allows Easy Transfers
- Car Style Seat Belts
- Solid Quick Release PU Tyres

MURRAYSIUK www.murraysmedical.co.uk **aktiv** wheelchairs

John Cousins
Midlands & North
07810 555885

Alex Lees
Midlands & South
07587 007699

Ph: 028 918 77951
Email: sales@murraysmedical.co.uk
www.murraysmedical.co.uk

MURRAYSIUK
www.murraysmedical.co.uk

aktiv wheelchairs



X2 Lite

X6 Kids



X7 Tilt

X2 Lite

Available in seat widths - 16", 18" & 20" / 18st User Capacity.



STANDARD FEATURES:

- Lightweight Aluminium folding frame with cable brakes
- Half-folding backrest
- Solid transit 12" or 24" self propel soft solid wheels
- Convertible from transit to self propel
- Flip back/removable desk arms as standard
- Swing in/out legrests with flip up footplates
- Padded, wipe clean, breathable black upholstery
- 7" front castors
- Generous push button release safety belt
- Attendant brakes

X5 HD Heavy Duty

Available in seat widths - 22", 24" & 26" / 25st User Capacity



STANDARD FEATURES:

- Lightweight Aluminium folding frame
- Half-folding backrest
- Solid transit 12" or 24" self propel soft solid wheels
- Convertible from transit to self propel
- Multiple axle heights and settings incl. amputee
- Flip back/removable desk arms as standard
- Swing in/out legrests with angle adjustable footplates
- One piece clip on footboard
- Padded, wipe clean, breathable tension adjustable black upholstery & 2" Seat cushion
- 8" robust front castors
- Height adjustable anti tippers
- Generous push button release safety belt

OPTIONS:

- Elevating legrests
- Stump support
- Amputee settings
- Short Desk, Loop or Height Adjustable armrests

X6 Kids

Available in seat widths - 12" & 14" / 15.5st User Capacity



STANDARD FEATURES:

- Folding lightweight aluminium frame
- Height adjustable push handles
- Flip back, detachable, height adjustable armrests
- Swing away leg rest, calf strap and clip-on one piece footboard
- Black seat sling with breathable, spacer fabric cover
- Tension adjustable backrest with breathable, spacer fabric cover
- 6" front castors
- Self propel 22" soft solid wheels
- Height adjustable anti-tippers
- Push button seat belt

OPTIONS:

- Elevating legrests

X3 Pro

Available in seat widths - 16", 18" & 20" / 21.4st User Capacity.



STANDARD FEATURES:

- Lightweight Aluminium folding frame
- Half-folding backrest
- Flip back and height adjustable armrests
- Solid transit 12" or 24" self propel soft solid wheels
- Convertible from transit to self propel
- Multiple axle heights and settings incl. amputee
- Rear seat height adjustable
- Swing in/out legrests with angle adjustable footplates
- 7" front castors with settable castor housing
- Height adjustable anti tippers
- Generous push button release safety belt

OPTIONS:

- Elevating legrests
- Stump support
- Amputee settings

X5 Modular

Available in seat widths - 16", 18" & 20" / 22st User Capacity



- Lightweight Aluminium folding frame
- Half-folding, angle adjustable backrest
- Solid transit 12" or 24" soft solid self propel wheels
- Convertible from transit to self propel
- Multiple axle heights and settings incl. amputee and sport
- Front and rear seat height adjustable
- Flip back/removable desk arms as standard
- Swing in/out legrests with angle adjustable footplates
- Padded, wipe clean, breathable tension adjustable black upholstery & 2" Seat cushion
- Height adjustable 7" front castors with settable castor housing
- Height adjustable anti tippers
- Generous push button release safety belt

OPTIONS:

- Elevating legrests
- Stump support
- Amputee settings
- Short Desk, Loop or Height Adjustable armrests

X7 Comfort Tilt & Recline

Available in seat widths - 16-18", 18-20" & 20-22" / 23.5st User Capacity



STANDARD FEATURES:

- Lightweight Aluminium rigid frame with tilt in space & recline functions
- Available in 16" quick release soft solid wheels
- Width and depth adjustable
- Attendant and occupant brakes
- Height adjustable backrest and push handle
- Removable, height adjustable armrests
- Elevating & length adjustable leg rests supplied as standard
- Multi-adjustable headrest and lateral support
- Adjustable Headrest